


LAMP FOR THE DETECTION OF ISKNV

PI2013702063



NEED

- To control the spread of ISKNV
- To test fish before moving/exporting
- A reliable & easy to use diagnostic kit
- Clear & easy to read results

INVENTION

- Rapid & accurate ISKNV diagnostic kit
- High specificity
- Ready-to-use reaction mixture
- Clear qualitative result
- Single post-reaction step

BENEFITS


- User friendly - Can be easily used at fish farm
- Economical - No sophisticated equipment needed
- Time saving - Detects ISKNV within 1 hour

USEFULNESS & APPLICATION

- Detection of ISKNV in fish
- Fast & efficient ISKNV diagnostic kit
- Can be used for field surveillance

POTENTIAL CONSUMERS

- Fish quarantine centre
- Aquaculture farms
- Private companies handling diagnostic kits/services



Project Leader : Prof. Dato' Dr. Mohamed Shariff Mohamed Din
 Co. Researchers : Kuttichantran Subramaniam, Prof. Dr. Abdul Rahman Omar, Prof. Dr. Mohd Hair-Bejo and Assoc. Prof. Dr. Ong Bee Lee
 Faculty : Faculty of Veterinary Medicine
 Email : shariff@upm.edu.my / pshariff@gmail.com
 Tel : +603.8609.3917
 Expertise : Fish health